



Compensation, Negative Controls and the Optimized Choice of Reagents for Multicolor Flow Cytometry

Overview


- Compensation
 - Introduction
 - Prevention of compensation related artifacts
- Negative Controls
 - Fluorescence Minus One (FMO)
 - Isotype controls
 - Transfection controls
- Combination of Reagents for Multicolor Flow Experiments
- Characteristics of Fluorochromes

Overview

- Compensation
 - Introduction
 - Prevention of compensation related artifacts
- Negative Controls
 - Fluorescence Minus One (FMO)
 - Isotype controls
 - Transfection controls
- Combination of Reagents for Multicolor Flow Experiments
- Characteristics of Fluorochromes


Compensation: Introduction

- Photons reaching a specific detector are coming from
 - the fluorochrome specific for this detector
 - optical background (especially cell-type specific auto-fluorescence)
 - photon spill over from all fluorochromes present in experiment




Compensation: Introduction

- Compensation is a procedure to subtracts from all photons that reach a detector the
 - 1) nonspecific electronic signal from cellular autofluorescence
Application of “Instrument Settings” for one specific (type of) sample
 - 2) nonspecific electronic signal from fluorochrome spill over



Compensation: Introduction

- Compensation is a procedure to subtracts from all photons that reach a detector the
 - 1) nonspecific electronic signal from cellular autofluorescence
Application of “Instrument Settings” for one specific (type of) sample
 - 2) nonspecific electronic signal from fluorochrome spill over





Compensation: Introduction

- 1) Generate sample-specific “Application Settings” according to SD_{EN} given in CS&T Baseline Report
 - Refer to Handout #3: Generation of Application Settings
 - Refer to the Presentation: QC for Digital Instruments: BD CS&T™
- 2) Perform an Automated Compensation using the BD FACS Diva™ software
 - Refer to Handout #4: Automated Compensation




Compensation: Prevention of compensation related artifacts

- To avoid errors in compensation
 - 1) never change the PMT V of fluorescent parameters (Instrument Settings) after a compensation has been calculated
 - 2) only the “median-related” = “statistic-based” compensation will result in accurate compensation values: Never compensate “by eye”
 - 3) use only compensation controls that show in minimum equally high signal intensities than the highest signal you expect in your sample
 - 4) BUT, the signal intensities of the compensation controls have to lie within the maximum linearity range of the scales (refer to the CS&T baseline report)
 - 5) it is not recommended to re-use compensations if a tandem dye is included in the panel

Compensation:

Prevention of compensation related artifacts

- Digital Compensation is a software-based compensation in that linear numbers are subtracted



Uncompensated List-Mode Data

	Time	FSC	SSC	FITC	PE
Event 1	0	60	120	5	23
Event 2	10	160	65	89	675
Event 3	30	650	160	39.271	30.621

- But if each optimal compensation is a subtraction that results in “0”, why is a compensation value of 20% worse than 0.2%?



Compensation:

Prevention of compensation related artifacts


- High compensation values induce “data-spread” effects that reduce substantially the resolution between populations

Uncompensated

Compensated



Gated on CD45 – FITC



Gated on CD45 – APC


CD3 – APC-H7

Overview

- Compensation
 - Introduction
 - Prevention of compensation related artifacts
- Negative Controls
 - Fluorescence Minus One (FMO)
 - Isotype controls
 - Transfection controls
- Combination of Reagents for Multicolor Flow Experiments
- Characteristics of Fluorochromes


Negative Controls: Fluorescence Minus One (FMO)

- Indicates the data-spread induced by “high” compensation values



Negative Controls: Fluorescence Minus One (FMO)

- Fluorochrome-combinations with “high” compensation values?




Fluorochrome	- % Fluorochrome	Spectral Overlap
PE	FITC	13.59
PerCP-Cy5-5	FITC	1.90
PE-Cy7	FITC	0.21
APC	FITC	0.01
APC-H7	FITC	0.00
horizon V450	FITC	0.00
horizon V500	FITC	0.59
FITC	PE	1.51
PerCP-Cy5-5	PE	16.44
PE-Cy7	PE	1.78
APC	PE	0.03
APC-H7	PE	0.00
horizon V450	PE	0.25
horizon V500	PE	0.21
FITC	PerCP-Cy5-5	0.00
PE	PerCP-Cy5-5	0.00
PE-Cy7	PerCP-Cy5-5	20.84
APC	PerCP-Cy5-5	0.34

Fluorochrome	- % Fluorochrome	Spectral Overlap
PE	FITC	25.64
PerCP-Cy5-5	FITC	2.39
PE-Cy7	FITC	0.49
APC	FITC	0.00
APC-Cy7	FITC	0.00
V450	FITC	0.50
V500	FITC	2.70
FITC	PE	0.82
PerCP-Cy5-5	PE	11.53
PE-Cy7	PE	2.00
APC	PE	0.02
APC-Cy7	PE	0.00
V450	PE	0.00
V500	PE	0.00
FITC	PerCP-Cy5-5	0.00
PE	PerCP-Cy5-5	0.00
PE-Cy7	PerCP-Cy5-5	37.06
APC	PerCP-Cy5-5	2.09


- Get experience: compare of multicolor compensation matrices.

Negative Controls: Fluorescence Minus One (FMO)

- Does every fluorochrome combination with “high” compensation values requests an FMO Control?




- Not if the fluorochromes bind to different cells!




Negative Controls: Isotype Control

- Indicates specificity of antigen-expression in unknown samples
- CD4 is THE T-helper cell marker**





- Indicates the nonspecific binding of antibodies to membranes





0.025 μ g antibody per 100 μ l whole blood

Negative Controls: Isotype Control

- Indicates specific binding of antibodies to very low and/or “smeared” expressed antigens



- But (if available and known) “biological controls” are better




Negative Controls: Isotype Control

- Isotype controls **HAVE** to
 - show perfect match in isotype subtype
 - come from the same species
 - be labeled to the same fluorochrome
 - show equalized concentrations between specific antibody and isotype control
 - be specifically purified from free fluorochromes for intra-cellular staining
 - show comparable antibody to fluorochrome ratios

Negative Controls: Transfection Control (“Mock Control”)

- Transfection methods using Lipofectin and similar substances enhance the autofluorescence due to the binding of plasmid-DNA to cell membranes
- Nontransfected cells are for this reason never a proper control!



- In this example the transfection-efficiency is 0% as the plasmid-promotor shows a deletion that prevents gene-expression!

Overview

- Compensation
 - Introduction
 - Prevention of compensation related artifacts
- Negative Controls
 - Fluorescence Minus One (FMO)
 - Isotype controls
 - Transfection controls
- Combination of Reagents for Multicolor Flow Experiments
- Characteristics of Fluorochromes



Combination of Reagents for Multicolor Flow Experiments

- Basic rule for the selection of antibody-fluorochrome-combinations:

Match the fluorochrome brightness with the antigen-expression levels


- For further information see in Appendix 4:

http://www.bdbiosciences.com/documents/Multicolor_AppNote.pdf

Combination of Reagents for Multicolor Flow Experiments

- Match fluorochrome brightness with antigen-expression levels
 - Fluorochrome brightness?

- The “Stain Index”



$$\text{Stain Index (SI)} = \frac{D}{W}$$

D = Δ between positive and negative peak medians

W = Spread of the background peak

- Stain indices of BD specific fluorochromes obtained on a BD FACSCanto™ II

Reagent	Clone	Stain Index
PE	RPA-T4	356.3
Alexa 647	RPA-T4	313.1
APC	RPA-T4	279.2
PE-Cy7	RPA-T4	278.5
PE-Cy5	RPA-T4	222.1
PerCP-Cy5.5	Leu-3a	92.7
BD Horizon™ V450	RPA-T4	90.0
Alexa 488	RPA-T4	75.4
FITC	RPA-T4	68.9
PerCP	Leu-3a	64.4
APC-Cy7	RPA-T4	42.2
Alexa 700	RPA-T4	39.9
AmCyan	RPA-T4	24.2
BD™ APC-H7	RPA-T4	18.0
BD Horizon™ V500	RPA-T4	12.0


Bright fluorochromes

Intermediate-bright fluorochromes

Dim fluorochromes

Combination of Reagents for Multicolor Flow Experiments

- Match fluorochrome brightness with antigen-levels
 - Antigen Expression level?





Combination of Reagents for Multicolor Flow Experiments

- Match fluorochrome brightness with antigen-levels
 - Antigen Expression level?

Or use literature-data about # of antigen-molecules per cell

Antigen	Antigen-Density	Expression Level
CD3	80.000	++
CD4	100.000	+++
CD14	144.000	+++
CD19	18.000	++
CD25	3.000	+
CD45	200.000	+++
CD56	10.000	+
CD127	2.000	+

Antigen-expression High / Intermediate / Low:

+++ / ++ / +

Example for one proper Antibody-Fluorochrome-Match for a BD FACS Aria™ 405nm / 688nm / 635nm

BD Horizon™ V450

FITC

BD™ APC-H7

PerCP-Cy™5.5

APC

BD Horizon™ V500

PE-Cy™7

PE

Overview

- Compensation
 - Introduction
 - Prevention of compensation related artifacts
- Negative Controls
 - Fluorescence Minus One (FMO)
 - Isotype controls
 - Transfection controls
- Combination of Reagents for Multicolor Flow Experiments
- Characteristics of Fluorochromes

Characteristics of Fluorochromes


- Select fluorochromes according to instrument configuration

Max Emission	Possible fluorochrome-excitation by			
	Violet 405nm	Blue 488nm	Y-G 561nm	Red 635nm
448nm or 452nm	BD Horizon™ V450 or Pacific Blue™			
461nm				
491nm	AmCyan or BD Horizon™ V500			
519nm or 520nm		AF®488 or FITC		
578nm		PE	PE	
660nm or 668nm				APC or AF®647
667nm or 678nm or 695nm		Pe-Cy™5 or PerCP or PerCP-Cy™5.5	Pe-Cy™5 or PerCP or PerCP-Cy™5.5	
723nm				AF®700
785nm or 783nm				APC-Cy™7 or BD™ APC-H7
785nm		Pe-Cy™7	Pe-Cy™7	



Characteristics of Fluorochromes


- Select fluorescent dyes according to instrument configuration



Max Emission	Possible fluorochrome-excitation by			
	Violet 405nm	Blue 488nm	Y-G 561nm	Red 635nm
461nm or 470nm or 477nm	Dapi or SYTOX-Blue or ECFP			
505 or 510nm		CFSE or EGFP		
522nm or 524nm or 528nm		SYBR-Green or SYTOX-Green or EYFP		
562nm			mOrange	
581nm			mTomato	
586nm		dsRed		
610nm			mCherry	
617nm or 647nm		PI or 7-AAD		
650nm			mPlum	
658nm				SYTOX Red
683nm				DRAQ5

Characteristics of Fluorochromes

- BD-supported fluorochromes





Characteristics of Fluorochromes: The “blue laser fluorochromes (I)”

	Advantages	Disadvantages
1. FITC	<ul style="list-style-type: none">▪ Most widely used▪ Stable, long-lasting conjugates▪ Easy conjugation (3-5 per Ig)▪ Cheapest fluorochrome	<ul style="list-style-type: none">▪ One of the dullest fluorochromes▪ Photobleaching (relevant only for microscopy)▪ pH sensitive
2. Alexa Fluor® 488	<ul style="list-style-type: none">▪ Very photo- and pH-stable (superior for microscopy)	<ul style="list-style-type: none">▪ One of the dullest fluorochromes
3. PE	<ul style="list-style-type: none">▪ Bright reagents: IC/FCM▪ Low backgrounds: IC/FCM 1st choice for intracellular protein/cytokine detection!▪ Good for quantification: 1 PE per Ig (BD™QuantiBrite)	<ul style="list-style-type: none">▪ Large size (240.000 D) is NOT a disadvantage!

Characteristics of Fluorochromes: The “blue laser fluorochromes (III)”

	Advantages	Disadvantages
6. PE-CyTM5	<ul style="list-style-type: none"> ▪ Bright fluorochrome ▪ Useful additional fluorochrome for the LSR II 	<ul style="list-style-type: none"> ▪ Unspecific binding to Fc Receptors ▪ Very strong spill over to APC: has to be adjusted often by titration ▪ Lot to lot differences in compensation ▪ Light sensitive (just keep samples in the dark)
7. PE-CyTM7	<ul style="list-style-type: none"> ▪ Bright fluorochrome: Preferable above PE-Cy5 	<ul style="list-style-type: none"> ▪ Lot to lot differences in compensation ▪ Light sensitive (just keep samples in the dark)



Characteristics of Fluorochromes: The “blue laser fluorochromes (II)”

	Advantages	Disadvantages
4. PerCP	<ul style="list-style-type: none">▪ Minimal Spill over to PE	<ul style="list-style-type: none">▪ Dull (equivalent to FITC)▪ Sensitive to photo-bleaching (relevant only for microscopy or high energy blue lasers. See BD FACSTTM Vantage)
5. PerCP-CyTM5.5	<ul style="list-style-type: none">▪ Intermediate bright fluorochrome (PerCP is dim)	<ul style="list-style-type: none">▪ Lot to lot differences in compensation▪ Light sensitive (just keep samples in the dark)



Characteristics of Fluorochromes: The “red laser fluorochromes (I)”

	Advantages	Disadvantages
1. APC	<ul style="list-style-type: none">▪ Bright fluorochrome▪ Low background: 2nd choice for intracellular protein / cytokine detection!	<ul style="list-style-type: none">▪ Large size (106.000 D) is NOT a disadvantage!
2. Alexa Fluor® 647	<ul style="list-style-type: none">▪ Equivalent to APC, BUT▪ Very photo- and pH-stable (superior for microscopy)	/
3. Alexa Fluor® 700	<ul style="list-style-type: none">▪ Very photostable (superior for microscopy)▪ Additional fluorochrome for the red laser	<ul style="list-style-type: none">▪ Filter set not included in standard BD FACSAria™ configuration, but easy to obtain additionally

Characteristics of Fluorochromes: The “red laser fluorochromes (II)”

	Advantages	Disadvantages
4. APC-CyTM7	<ul style="list-style-type: none">▪ Additional red laser excited fluorochrome	<ul style="list-style-type: none">▪ Dim▪ Photo-instable▪ Fixative-sensitive▪ Lot to lot differences in compensation
5. BDTM APC-H7	<ul style="list-style-type: none">▪ Photo-stable▪ Fixative-stable: Preferable above APC-Cy7	<ul style="list-style-type: none">▪ Dim▪ Lot to lot differences in compensation



Characteristics of Fluorochromes: The “violet laser fluorochromes (I)”

	Advantages	Disadvantages
1. Pacific Blue™	<ul style="list-style-type: none">▪ Little Spill over to FITC	<ul style="list-style-type: none">▪ dim
2. BD Horizon™ V450	<ul style="list-style-type: none">▪ Little Spill over to FITC▪ Intermediate bright	/



Characteristics of Fluorochromes: The “violet laser fluorochromes (II)”

	Advantages	Disadvantages
3. AmCyan	<ul style="list-style-type: none">▪ Before BD Horizon™ V500 and BD Horizon™ V450 the only violet-excited fluorochrome from BD	<ul style="list-style-type: none">▪ dim▪ extremely strong spill over to FITC▪ Fixative-sensitive▪ Light-sensitive
4. BD Horizon V500	<ul style="list-style-type: none">▪ Little spill over to FITC: highly preferable alternative to AmCyan▪ Fixative-stable	<ul style="list-style-type: none">▪ Dim

Characteristics of Fluorochromes

- Summary:
 - There are no “bad” fluorochromes!
 - Considering limitations, all fluorochromes are useful designing multicolor flow experiments
 - But some fluorochromes are easier to handle for multicolor-combinations than others!