

BM Staining

- Isolate bone marrow from Bl6 mice:

1.) collect all legs needed for the experiment

2.) flush the bones with ice-cold MACS buffer into a 6 Well Plate (filled with 2ml MACS buffer; ice-cold)

- Filter Cells with a 70µm cell strainer; add additional MACS buffer and wash filter
- Centrifuge **(1500 rpm, 4°C, 10 min)**
- *Do Ery-Lysis only if necessary*
- *As some of the antibody-stainings are Ca²⁺ dependent EDTA could disturb your staining.*

If your staining panel works with EDTA, MACS buffer is recommended from this step on.

- Discard supernatant and resuspend in 1 ml **MACS/FACS** buffer (+ - **EDTA**)
- count cells
- concentrate cells to a maximum of 50Mio/ml in **MACS/FACS** buffer (+ - **EDTA**)
- Stain with primary antibodies for 45min at 4°C.
- Centrifuge and add **MACS buffer** (with EDTA) 1ml/30Mio Cells

Material:

Ery-Lysis-Buffer

0,15 M NH₄Cl
0,02 M HEPES
0,1 mM EDTA

MACS Buffer

1x PBS
2 mM EDTA
2% FCS

FACS Buffer

1x PBS
2% FCS

Cell strainer 70µm Nylon (BD Falcon #352350)

