

## Cell isolation protocol for T reg sorting

(with/out enrichment)

Protocol based on stroma isolation protocol: Fletcher, AL et al. (2011), *Frontiers in Immunology*

### Required reagents:

- 1x PBS, pH 7.4
- 1x PBS SE, pH 7.4 (PBS, 2%(v/v) FCS, 1-5mM EDTA (as you wish))
- RPMI, 2%(v/v) FCS
- Collagenase D or P (both work equally, from Roche 11088882001)
- DNase I, grade II (from Roche 10104159001)
- 0.5M EDTA
- PBC lysis buffer (Biolgend 420301) or any ACK Buffer

### Procedure:

- Prepare fresh digestion Medium:

RPMI, 2%(v/v) FCS  
+ 0.2mg/ml Collagenase D  
+ 0.1mg/ml DNase I

15ml/mouse, keep at 37°C in a water bath

- Sacrifice mice, take out spleen/LNs, put into ice cold PBS (**without** EDTA!)
- Put spleen/LNs into a petri dish. Cut spleen/LNs into 4-6/2-4 pieces and transfer into a fresh 15ml falcon  
*(To make things easier: if you are only interested in LNs, a 1.5ml reaction tube is sufficient. Incubate then at 1000rpm on a thermo-shaker, 37°C)*
- Add 3ml of 37°C Digestion buffer & incubate for 10min, shake tube every 3min

- After 10min, medium should change from clear to turbid:  
Then put the **Supernatant** (containing the already dissociated cells) into a **fresh 50ml falcon**, add to the SN EDTA qs 5mM (e.g. for 2.5ml, 25 $\mu$ l of 0.5M EDTA), mix and spin 5min at 380g. Discard the digestion medium and resuspend in normal and ice cold RPMI 2%(v/v) FCS, keep on ice.  
At the same time, refill **the digestion tube** with 3ml of digestion buffer and put it back into the 37°C water bath, shake tube every 3min.
- After another 10min, try to pipette smoothly the debris with a p1000 (cut the tip)  
Take out the SN and pool it to the already isolated cells. Add EDTA qs 5mM and spin and resolve pellet in ice cold RPMI 2%(v/v) FCS. Keep on ice.  
At the same time, add 3ml of fresh digestion buffer to the digestion tube containing the organ debris.  
  
  - ⊆ Repeat this **3 more times** (total 60-70min). All debris should be digested after 50-70min. Do not forget to pipette every 5min with a cut tip.
- Adjust the volume to 50ml RPMI 2%(v/v) FCS, filter isolated cells through a 70 $\mu$ m cell strainer, spin 10min at 380g (washing step/filtering step)
- Resuspend the pellet in 5ml of 1x RBC lysis buffer, incubate for 4-5min, add PBS SE qs 30ml, spin at 380g, 10min.
- Discard supernatant, pellet should be brownish/white. Resolve pellet in 10ml per mouse in PBS SE. Eventually filter again.
- Count cells in Trypan blue (pre-dilute 1:100)
- Adjust volume to 1ml per 20 x 10<sup>6</sup> cells with PBS SE
- Add Fc-block (anti-CD16/CD32, 2.4G2, BD Pharmingen) at 1 $\mu$ g/ml (1/500) and incubate for 5min on a roller shaker at RT
- Without wash, proceed to staining:  
  - CD4-A647** (RM4-5, Biolegend) at 0.33 $\mu$ g/ml (1/1500)
  - CD25-PE** at 0.4 $\mu$ g/ml (PC61, Biolegend; not for in vivo, depleting via ADCC)  
 Incubate 10-15min on a roller shaker at 4°C
- Spin at 380g, 10min, discard SN, fill tube with PBS SE, spin at 380g, 10min, discard SN
- Take up cells in 40 x 10<sup>6</sup> cells per ml
- Proceed to sort

Scheme-Digestion:

