

Pre-pro-B Sort

- Isolate bone marrow from Bl6 mice
 - Resuspend cells and centrifuge (**1500 rpm, 4°C, 1min per 2ml**)
 - Discard supernatant
 - Resuspend in Ery-Lysis-buffer (2ml/mouse) for 5min at RT
 - Stop lysis by adding MACS buffer (2ml/mouse) and centrifuge
 - Discard supernatant
 - Resuspend in 5 ml MACS buffer and filter (70µm); add 5ml MACS buffer and wash filter; take aliquot and count cells.
 - Centrifuge and discard supernatant
 - concentrate cells to 10Mio/ml in MACS buffer
 - FC block (1:100)
 - Centrifuge and discard supernatant
 - concentrate cells to 10Mio/ml in **FACS buffer** (without **EDTA**)
 - Split Cells:
- (Titration: 6 x antibodies x 5 titration steps: 30 vials, each 200.000 cells; total: 6Mio)
- Single staining: 6x 200.000 Cells; total 1.2 Mio; each in 200 µl vol.
 - Take remaining cells for panel staining
 - Stain with primary and Biotin antibodies for 30min at 4°C.
 - Centrifuge and add **FACS buffer** 1ml/10Mio Cells
 - Stain with secondary antibody (SA) for 30min at 4°C.
 - Centrifuge and add MACS buffer 1ml/10Mio Cells

Material:

Ery-Lysis-Buffer

0,15 M NH₄Cl
0,02 M HEPES
0,1 mM EDTA

MACS Buffer

1x PBS
2 mM EDTA
2% FCS

FACS Buffer

1x PBS
2% FCS

Cell strainer 70µm Nylon (BD Falcon #352350)

Panel staining:

Primary staining:

AA4.1 (CD93)-FITC
cKIT-PE
B220-PerCP
HSA (CD24)-PE-Cy7
CD43-APC
CD19-APC-Cy7
Cd11b Biotin (dump channel)
CD3 Biotin (dump channel)
Ly6G Biotin (dump channel)

Secondary staining:

Streptavidin A350

Staining of prepro B Cells:

AA4.1 POS
CKIT POS
B220 POS
HSA Medium
CD43 POS
CD19 NEG
Dump channel NEG